


# Reprocessing quick reference

## Elements of a designated reprocessing area



## Elements of reprocessing reusable instruments

**Personal Protective Equipment (PPE) must always be available**

Masks, gowns, gloves (puncture-resistant) and eye protection are necessary when splashes or sprays may occur.

| Pre-cleaning                                    | Transport                                                                                                                          | Disassembly                                                                                                                         | Cleaning                                                                                                                                                                                           | Drying                                                                                            | Sterilization                                                                                                                                                                                                                                                                                                                                              |
|-------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Remove gross soil immediately to prevent drying | Transport container must be cleanable, closed, covered and puncture proof<br>Dispose of sharps promptly (ideally at point of care) | Disassemble as per Manufacturer's Instructions for Use (MIFU)<br>Inspect for damage<br>Sort instruments into sets<br>Soak/pre-treat | Completely submerge equipment<br>Manual friction with detergent or mechanical cleaning if available (e.g. ultrasonic cleaner or washer/disinfector)<br>*mechanical cleaning is preferred.<br>Rinse | Use lint-free towel<br>Compressed filtered air to dry lumens<br>Lubricate, reassemble as required | Items must not touch in package, open hinged items<br>Wrap/label with date, sterilizer used, load number, and initials. Add load contents if not visible.<br>Ensure quality indicators are documented: <ul style="list-style-type: none"><li>Physical parameters</li><li>Chemical indicator – internal and external</li><li>Biological indicator</li></ul> |

### Quality control

#### Sterilization process

- ✓ Physical indicators (time, temperature, pressure) – monitor with each load by the sterilizer printout or USB
- ✓ Chemical indicators – monitor on each package
  - Internal indicators should be a minimum Class
  - External indicators are Class 1
  - For dynamic air removal type sterilizer, perform a Class 2 indicator (Bowie-Dick) once daily in a process challenge device (PCD)
- ✓ Biological indicators – monitor once daily minimum
  - Place a biological indicator in a PCD in the first load of each sterilizer in use daily
  - Use a biological indicator in every load containing implants
  - Instruments must be held until the result of the biological indicator is available. To use the instruments prior to results, an internal Class 5 or 6 chemical indicator should be in a PCD in each load

#### Ultrasonic cleaners

- ✓ Performance testing done weekly according to manufacturer's instructions

#### High level disinfectant

- ✓ Use test strips each day it is used according to manufacturer's instructions

**DOCUMENTATION REQUIRED**

**Single-use medical/dental instruments and equipment are usually labelled by the manufacturer with a symbol.**



Any equipment labelled as single-use can never be reprocessed or re-used. Examples: syringes, needles, some suction tips, etc.

## Storage

**Never store sterile/clean instruments in the soiled area of the reprocessing room, on the floor, under/adjacent to sinks or with unclean items.**

#### Store single-use disposable semi-critical items

- Store in a clean, dry, covered area and handled with clean hands or forceps

#### Store reprocessed sterile or critical items

- In sterile package until time of use
- in clean, dry, secure area that prevents contamination
- Only handle with clean hands
- Storage is not accessible from soil reprocessing area
- Ensure items that have been reprocessed can be differentiated from items that have not (i.e. colour coding)

## Key principles for reprocessing instruments

It is essential to maintain an overall inventory of instruments used in each facility.

### The inventory must have:

- ✓ The recommended reprocessing instructions for each instrument
- ✓ Reflect the manufacturer's instructions for use for each item
- ✓ Be shared and accessible to all staff involved in reprocessing

## Policies need to include:

### Cleaning and maintenance

- Environmental cleaning process and schedule for reprocessing area
- Scheduled maintenance for cleaning and sterilization equipment

### Reprocessing

- All aspects of reprocessing
- Process for recall of instruments if there is a reprocessing failure
- Quality monitoring and documentation
- Items are reprocessed as per manufacturer's instructions or are single use
- Single-use items are not reprocessed

### Occupational Health and Safety


- Response to blood/body fluid spills or exposures
- Prevention/management of sharps injuries
- Use of personal protective equipment
- All chemicals are approved by Health Canada

### Education and training

- Specify requirements for staff education and training
- Frequency of education and training at regular intervals (e.g. orientation and continuing education)
- Competency assessment for all personnel involved in reprocessing



## Reprocessing decision tree



MIFU = Manufacturer's Instructions for Use



If there is a discrepancy between the reprocessing level recommended by the manufacturer and the intended use of the instrument by Spaulding's criteria, the higher level of disinfection/sterilization must be used.\*\*

## Helpful resources

Public Health Ontario – [Reprocessing resources](#)

Public Health Ontario – [IPAC online learning](#)

Public Health Ontario – [Resources for preventing IPAC lapses](#)

[Best Practices for Cleaning, Disinfection and Sterilization of Medical Equipment/Devices in all Health Care Settings](#) (includes Spaulding's Criteria)



\*\*If the MIFUs not compatible with Spaulding's criteria, alternative device or instrument should be considered.

Adapted and reprinted with permission of the Region of Waterloo Public Health and Paramedic Services.